



MULTINFOR

Equipamento e Manutenção Hospitalar Lda.

Q5/Q3

Modular Patient Monitor



NEW

TECHNICAL SPECIFICATION:_____Q5/Q3

Q5/Q3 Main Unit

Size and Weight

- Size : 318mmX264mmX152mm
- Weight : ≤4.5kg
- Standard module slot: 2

- Power supply
- Power Voltage: AC 100-240V 50/60Hz
 - Input Current: 1.1~0.5A
 - Safety class: Category I

- Display
- Q5 : 12.1"/Q3:10.4" Color TFT-LCD
 - Resolution : 800X600 pixels

- Battery
- Type: Rechargeable Lithium battery, 11.1V/4.0AH
 - Operating time under the normal use and full charge : ≥210minutes

- Recorder (Option)
- Method : Thermal dot array
 - Paper width : 50 mm (1.97 in)
 - Paper length: 15m
 - Paper Speed : 12.5 / 25 / 50 (mm/sec)
 - Traces : Maximum 3 tracks
 - Recording way : Real-time recording, Periodic recording, Alarm recording

- Alarm
- Level : Low, medium and high
 - Indication : Auditory and visual
 - Patient Physiological Alarm Light color: Yellow & Red:
 - Equipment Technical Alarm Light color: Blue
 - Supports Pitch Tone and multi-level volume;
 - Supports custom arrhythmia tone

- Input device
- Touch screen:standard config
 - Knob:standard config
 - Mouse input:Support
 - Keyboard input:Support

- System Output & Extensible Interfaces
- Ethernet Network: 1 Standard RJ45 socket
 - Defibrillation Output:1 BNC connector
 - Nurse Call: 1 RJ11 connector
 - Video Output : 1 VGA port
 - USB1.1 port : 2
 - SD meomery card : 2G (Option)
 - Analog Output (ECG or IBP) : Option

- Trend & Reviewing :
- Trend : Lond trend:168h,minimum resolution is 1min (store when power goes off)
High resolution trend:2h, minimum resolution is 5s
 - NIBP measurement reviewing : 1000 groups
 - ARR event: 128 groups of ARR event and the associated waveform.
 - Alarm events:128 groups of parameter alarm events and associated parameter waveform at the alarm moment
 - Full Disclosure waveform: 24 hours for 3 waveforms(with 2G SD cord)

- Environment
- Operating temperature:0℃ ~ +40℃
 - Storage temperature:-20℃ ~ +50℃
 - Operating humidity:15% to 85% non condensing
 - Storage humidity:10% to 93% non condensing
 - Operating atmospheric pressure: 860hPa to 1060hPa
 - Storage atmospheric pressure:500hPa to 1060hPa

- Safety:
- IEC60601-1 Approved, CE marking according to MDD93/42/EEC

Performance:

- ECG
- Lead Mode : 3-leads ECG input
5-leads ECG input
12-leads ECG input
 - Lead selection :I, II, III
I, II, III, aVR, aVL, aVF, V-
I, II, III, aVR, aVL, aVF, V1~V6 (option)
2.5 mm/mV(x0.25), 5 mm/mV(x0.5), 10 mm/mV(x1),
20 mm/mV(X2),40mm/mv(x4),Auto
 - CMRR : Monitor mode≥105dB
Surgery mode ≥105dB
Diagnostic mode≥90dB
 - Frequency response (-3d B):
Monitor mode 0.5~40Hz
Surgery mode 1~25Hz
Diagnostic mode 0.05~150Hz
 - Input impedance : ≥5.0 Mohm
 - ECG signal range: ± 10.0mV
 - Electrode offset potential: ± 500mV
 - Patient Leakage Current : <10 uA
 - Standardizing signal :1 mV ± 5%
 - Baseline recovery : <5s after Defibrillation. (Mon or Surg mode)
 - Indication of electrode separation : Every electrode (exclusive of RL)
 - Protection: Breakdown Voltage 4000VAC 50/60Hz; defibrillator proof
 - Sweep speed : 12.5mm/s, 25mm/s, 50mm/s

- HR
- Range : Adult 10~300 bpm
Pediatric & Neonate: 10 ~350bpm

- Refreshing time : ≤50 bpm Per 2 pulses
50~120bpm Per 4 pulses
≥120bpm Per 6 pulses
- Resolution : 1 bpm
 - Accuracy : ± 1% or ± 1bpm, whichever is greater

- ST segment
- Measurement range : -2.0mV~2.0mV
 - Accuracy: -0.8mV~0.8mV ; ± 0.02mV or ± 10%, whichever is greater
Over±0.8mV: unspecified
 - Resolution : 0.01mV

- RESP
- Method : Thoracic impedance
 - Lead Selected from: I (RA-LA) or II (RA-LL); Default: I
 - Gain : x0.25, x1 x2 x4
 - Bandwidth: 0.25 Hz to 2.0Hz (-3dB)
 - Sweep speed : 6.25mm/s, 12.5mm/s, 25mm/s
 - Measurement Range : 0~150 rpm
 - Resolution : 1 rpm
 - Accuracy : ± 2rpm or ± 2%, whichever is greater
 - Delay of Apnea Alarm : 10s,15s,25s,30s,35s,40s,45s,50s,55s, 60s

- NIBP
- Way of measurement : Automatic oscillometry
 - Range of measurement :
Adult : SYS 30~270 mmHg
DIA 10~220 mmHg
MAP 20~235 mmHg
Child: SYS 30~235 mmHg
DIA 10~220 mmHg
MAP 20~225 mmHg
Neonate: SYS 30~135 mmHg
DIA 10~100 mmHg
MAP 20~125 mmHg
 - Cuff pressure range : 0~300 mmHg
 - Resolution : 1 mmHg
 - Pressure Accuracy : Static : ± 2% or ± 3mmHg , whichever is greater
Clinical : ±5 mmHg average errorr
standard deviation : ≤8 mmHg
 - Unit: mmHg, kPa
 - Measurement mode:Manual , Auto, STAT
 - Intervals for AUTO measurement time : 1,2,3,4,5,10,15,30,60,90 minutes;
2,4,8,12hours
 - STAT mode cycle time : Keep 5 minutes, at 5 seconds interval.
 - Overpressure Protection : Hardware and software double protections
 - Pulse rate range : 40 ~ 240 bpm

- BLT-SpO2 (Digital Technic)
- Measurement Range : 0~100%
 - Resolution: 1%
 - Accuracy: At 70~100%, ±2%
At 0~69%, unspecified

- PR
- Measurement Range : 25~255 bpm
 - Resolution : 1 bpm
 - Accuracy : ± 1% or ± 1 bpm, whichever is greater

- Nellcor-SpO2 (option)
- Measurement Range : 0~100%
 - Resolution : 1%
 - Accuracy : At 70~100%, ±2% (Adult)
At 70~100%, ±3% (Neonate)
At 70~100%, ±2% (Low Perfusion)
At 0~69%, unspecified

- PR
- Measurement Range : 20~300 bpm
 - Resolution : 1 bpm
 - Accuracy : 20bpm to 250bpm: ± 3 bpm
251bpm to 300bpm: unspecified

- Masimo SpO2 (option)
- Measurement range:0% to 100%
 - Resolution: 1%
 - Accuracy: 70% to 100% ± 2% Adult/pediatric, non-motion conditions
70% to 100% ±3% Neonate, non-motion conditions
70% to 100% ±3% Motion conditions
0% to 69% unspecified
 - Average time: 2-4s, 4-6s, 8s, 10s, 12s, 14s, 16s

- PR
- Measurement range:25 bpm to 240 bpm
 - Accuracy : ± 3bpm Non-motion conditions
± 5bpm motion conditions
 - Resolution:1 bpm

- TEMP
- Max Channel : 8
 - Measurement way: Thermal resistance way
 - Measurement Range : 0.0℃ ~ 50.0℃ (32°F~122°F)
 - Accuracy : ±0.1℃ or ±1 ° F (exclusive of probe)
 - Resolution : 0.1℃ or 1 °F
 - Unit : Celsius (℃), Fahrenheit (° F)

- IBP
- Max Channel : 8
 - Measurement way: Directly invasive pressure measurement
 - Sensitivity of transducer: 5uV/V/ mmHg, ± 2%
 - Impedance of transducer: 300 to 3000Ω

- Measurement Range : -50 ~ +350 mmHg
- Resolution : 1mmHg
- Unit : mmHg, kPa,cmH2O
- Accuracy :
Static: ± 1mmHg or 2%, whichever is greater (exclusive of transducer)
± 4mmHg or 4%, whichever is greater (inclusive of transducer)
Dynamic : ± 4mmHg or 4%, whichever is greater
Transducer sites : Arterial Pressure (ART)
Pulmonary Artery Pressure (PA)
Left Atrium Pressure (LAP)
Right Atrium Pressure (RAP)
Central Venous Pressure (CVP)
Intracranial Pressure (ICP)
P1/P2
- Selection of measurement range :
ART : 0 ~ +350mmHg
PA : -10 ~ +120 mmHg
CVP/RAP/LAP/ICP: -10 ~ +40 mmHg
P1/P2 : -50 ~ +350 mmHg

- EtCO₂ (Sidestream)
- Measure method : Infrared spectrum
 - Measurement Range : 0.0~13.1% (0~99.6 mmHg)
 - Resolution : 1 mmHg
 - Unit : %, mmHg, kPa
 - Accuracy : 0% to 4.9% ± 0.3% (± 2mmHg)
5.0% to 13.1%, < ± 10% of reading
 - Measurement range of awRR : 3~150 rpm
 - Calibration: Offset calibration: auto, manual, Gain calibration

- EtCO₂ (Mainstream)
- Measure method : Infrared spectrum
 - Warm up time : Capnogram displayed in less than 15 seconds, At an ambient temperature of 25℃ ,full specifications within 2 minutes.
 - Measurement Range : 0.0~19.7% (0~150 mmHg)
 - Resolution : 1 mmHg
 - Rise time (10 l/min) : ≤ 60 ms
 - Unit : %, mmHg, kPa
 - CO₂ Accuracy 0 - 40 mmHg, ± 2mmHg
41 - 70 mmHg, ± 5% of reading
71 - 100 mmHg, ± 8% or reading
101 - 150 mmHg, ± 10% of reading
(at 760 mmHg, ambient temperature of 35℃)
 - awRR measurement range: 0~150 rpm
 - awRR measurement Accuracy: ± 1 rpm

- EtCO₂ (Microstream)
- Measure method : Infrared spectrum
 - Warm up time : Capnogram displayed in less than 20 seconds, At an ambient temperature of 25℃ , full specifications within 2 minutes.
 - Measurement Range : 0 - 19.7%(0-150mmHg)
 - Resolution : 1mmHg
 - Unit : %, mmHg, kPa
 - CO₂ Accuracy 0 - 40 mmHg, ± 2mmHg
41 - 70 mmHg, ± 5% of reading
71 - 100 mmHg, ± 8% of reading
101 - 150 mmHg, ± 10% of reading
(at 760 mmHg, ambient temperature of 25℃)
(when Rr >80 rpm, all the range is ± 12% or reading)
CO₂ response time: < 3s
 - awRR measurement range 2~150 bpm
 - awRR measuremen Accuracy : ± 1rpm
 - Sample Flow Rate 50 ml/min ± 10ml/min

- Anesthetic Gas
- Measure method : Infrared spectrum
 - Measure mode : Mainstream or Sidestream
 - Fi and Et values : CO₂ N₂O O₂ AG (HAL, ISO, ENF, SEV, DES)
 - Resolution : 1%
 - Unit : %.
 - Calibration : Room air calibration performed automatically when changing airway adapter (<5 sec)
 - Warm-up time : <10 s, full accuracy within 1 min
 - Measurement and alarm range of AG:

Gas	Range	Accuracy
CO ₂	0-10 %	± (0.3%ABS+4% REL)
N ₂ O	0-100 %	± (2% ABS+8% REL)
O ₂	10-100 %	± (2% ABS+2% REL)
HAL, ISO, ENF	0-5%	± (0.15% ABS+10% REL)
SEV	0-8%	± (0.15% ABS+10% REL)
DES	0-18%	± (0.15% ABS+10% REL)
- awRR measurement range :	0~150 rpm	
- awRR measurement Accuracy :	± 1 rpm	
- Rise time (flowing speed 10 l/min)	CO ₂ ≤ 90 ms O ₂ ≤ 300 ms N ₂ O ≤ 300 ms Hal, Iso, Enf, Sev, Des ≤ 300 ms	
 - Total system response time : < 1 seconds

- C.O.
- Measurement Mode: Thermal dilution method
 - Measurement Wave: Thermal dilution curve
 - Measurement parameters: C.O., TB, TI, C.I.
 - Measurement Range:C.O.: 0.1 L/min ~ 20 L/min
TB: 23.0 ~ 43.0℃
TI: -1.0 ~ 27.0℃
 - Resolution :C.O.: 0.1 L/min
TB: 0.1℃
TI: 0.1℃
 - Accuracy:C.O.: 2% SD TB, TI: ±0.1℃
 - TB Alarm range:23.0~43.0℃ , high/low limit can be adjusted continuously

Standard configuration:

Mainunit: Q5 : 12.1"/Q3:10.4" TFT-LCD display, 2 Standard module slot, Touch Screen, 1 RJ45 enthernet socket, 1 Defibrillation Output, 1 Nurse Call socket, 1 VGA port, 2 USB1.1 port, 1 Lithium rechargeable battery.

Option:

Option Module: Sidestream CO₂ Module, Microstream CO₂ module, Mainstream CO₂ module, AG module, C.O. module, IBP module, Temp module, Masimo SpO₂ module, Nellcor SpO₂ module.
Navigating: USB compatible mouse and keyboard.
Printing: 3 channel thermal recorder
Mounting: Rolling stand , wall mount
Battery: 11.1V/4.0AH Rechargeable Lithium Battery.
Other options: External Display, Wireless Lan, Extensive Memory card, Analog Output (ECG or IBP)



*Specifications subject to change without prior notice.

BIOLIGHT CO.,LTD.

Add: NO.2 Innovation First Road, Technical Innovation Coast, Hi-tech Zone, Zhuhai, P.R.China Postcode: 519085

Tel: +86-756-3399935 3399999 Fax: +86-756-3399911 E-Mail: Overseas@blt.com.cn



Http: //www.blt.com.cn

Q5/Q3-Y1303-1